


## DMA No. 9

## Dynamic Viscoelastic Data of Polyethyleneterephthalate

**1. Sample** Polyethyleneterephthalate : PET  
 (Brand Name : S-Type film, biaxially oriented)

**2. Chemical Structure**



**3. Thermal History** Annealed at 130 for 1 hour.

**4. Instruments** SDM5600 Rheol. Station  
 DMS200 Dynamic Mechanical Spectrometer

**5. Conditions**  
 Deformation mode : Tensile mode  
 Sample Size : 10.00(ℓ) × 10.00(w) × 0.01(t)mm  
 Temperature Range : -150 ~ 220  
 Heating Rate : 1K/min  
 Atmosphere : N<sub>2</sub>  
 Frequency : 0.5,1,2,5,10Hz

**6. Transition temperature and activation energy based on tanδ**

	Transition Temperature ( )	ΔEa (kJ/mol)	Comments
α Transition	114 (1Hz)	398	Glass Transition
β Transition	-70 (1Hz)	76	Local mode relaxation


**7. Thermal Analysis Data**

Annealed Sample (1st run) T<sub>m</sub> : 257.7 , ΔH<sub>m</sub> : 61.0 J/g  
 Quenched Sample (2nd run) T<sub>g</sub> : 76.5 , ΔC<sub>p</sub> : 0.377 J/deg·g  
 T<sub>c</sub> : 134.0 , ΔH<sub>c</sub> : 37.9 J/g  
 T<sub>m</sub> : 257.5 , ΔH<sub>m</sub> : 51.3 J/g

DSC 10K/min

DMS

Name: Sample: PET film  
 Date: 91/03/27 21:21 Temp. mode: Ramp  
 Comment: 1°C/min Deform: Tension  
 1\*s: 10.000 mm \* 0.100 mm2  
 Frequency: 0.5 ~ 10 Hz



DMS

Name: Sample: PET film  
 Date: 91/03/27 21:21 Temp. mode: Ramp  
 Comment: 1°C/min Deform: Tension  
 1\*s: 10.000 mm \* 0.100 mm2  
 Frequency: 0.5 ~ 10 Hz

